

## Product Information

**Kit Name:** VetPCR™ PBFVDV Detection Kit

**Storage Recommendations:** -20°C

VetPCR™ PBFVDV Detection Kit is the direct detection of Psitticine Beak and Feather Disease on the basis of a genetic database, so it can diagnose very fast and accurately.

It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

The kit consists of the following components that are enough for amplification of genomic DNA:

1. VetPCR™ PBFVDV PCR Pre-mixture
2. DNase/RNase free water
3. PBFVDV PCR Positive control
4. PBFVDV PCR Positive control Pre-Mixture
5. Brig™ Molecular Weight marker
6. Mineral Oil Solution.

**QUALITY CONTROL** - This batch meets the following release criteria: Each components of the system has been tested and passed its pH, conductivity or RNase activity measurements according to individual quality control specifications.

**Bioingentech Genomic DNA Purification Kit** is a registered trademark of Bioingentech Corporation.

**Unless other stated, the performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.**

### **PRODUCTS USE LIMITATIONS, WARRANTY, DISCLAIMER.**

Bioingentech manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products.

Bioingentech products contain chemicals which may be harmful if misused. Due care should be exercised with all Bioingentech products to prevent direct human contact.

Each Bioingentech product is shipped with documentation stating specifications and other technical information. Bioingentech products are warranted to meet or exceed the stated specifications. Bioingentech's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Bioingentech makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION OR ANY OTHER MATTER WITH RESPECT TO BIOINGENTECH PRODUCTS. In no event shall Bioingentech be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence), or strict liability arising in connection with the sale or the failure of Bioingentech products to perform in accordance with the stated specifications.

Bioingentech Ltd.

Salas 350, Piso 2

Concepción

CHILE

Telephone (+56) 41-2790435

Website [www.bioingentech.com](http://www.bioingentech.com)

Info [info@bioingentech.com](mailto:info@bioingentech.com)

2004 – 2009 BioInGenTech corporation.  
All rights reserved.