

Product Information

Kit Name: VetPCR™ AIV Detection Kit

Storage Recommendations: -20°C

VetPCR™ AIV Detection Kit is the direct detection of Avian influenza virus on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 4 hours for detection. Therefore, it is a very fast accurate and reliable technique.

The kit consist of the following components that are enough for amplification of cDNA from Total RNA:

1. VetPCR™ AIV RT-PCR Pre-mixture
2. VetPCR™ AIV PCR Pre-mixture
3. BrigRT-PCR™ solution
4. Biotech™ Transcriptase solution
5. DNase/RNase-free water
6. AIV RT-PCR Positive control
7. AIV RT-PCR Positive control Pre-mixture
8. AIV PCR Positive control Pre-mixture
9. Brig™ Molecular Weight marker
10. Mineral Oil solution

QUALITY CONTROL - This batch meets the following release criteria: Each components of the system has been tested and passed its pH, conductivity or RNase activity measurements according to individual quality control specifications.

VetPCR™ AIV Detection Kit is a registered trademark of Bioingentech Corporation.

Unless other stated, the performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

PRODUCTS USE LIMITATIONS, WARRANTY, DISCLAIMER.

Bioingentech manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products.

Bioingentech products contain chemicals which may be harmful if misused. Due care should be exercised with all Bioingentech products to prevent direct human contact.

Each Bioingentech product is shipped with documentation stating specifications and other technical information. Bioingentech products are warranted to meet or exceed the stated specifications. Bioingentech's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Bioingentech makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE. MERCHANTABILITY, CONDITION OR ANY OTHER MATTER WITH RESPECT TO BIOINGENTECH PRODUCTS. In no event shall Bioingentech be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence), or strict liability arising in connection with the sale or the failure of Bioingentech products to perform in accordance with the stated specifications.

Bioingentech Ltd.

Salas 350, Piso 2

Concepción

CHILE

Telephone (+56) 41-2790435

Website www.bioingentech.com

Info info@bioingentech.com

2004 – 2009 BioInGenTech corporation.
All rights reserved.