

Product Information

Kit Name: VetPCR™ BAB.BOVI Detection Kit

Storage Recomendations: -20°C

VetPCR™ BAB.BOVI Detection Kit is the direct detection of Babesia bovis on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

PRODUCTS USE LIMITATIONS, WARRANTY, DISCLAIMER.

Bioingentech manufactures products for a number for intendend uses. Please refer to the product label for the intendend use statements for specific products.

Bioingentech products contain chemicals which may be harmful if misused. Due care should be exercised with all Bioingentech products to prevent direct human contact.

Each Bioingentech products is shipped with documentation stating specifications and other technical information. Bioingentech products are warranted to meet or exceed the stated specifications. Bioingentech's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Bioingentech makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING WITHOUTH LIMITATION AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE. MERCHANTABILITY, CONDITION OR ANY OTHER MATTER WITH RESPECT TO BIOINGENTECH PRODUCTS. In no event shall Bioingentech be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whetever based upon warranty, contract, tort (including negligence), or strict liability arising in connection with the sale or the failure of Bioingentech products to perform in accordance whit the stated specifications.

The kit consists of the following components that are enough for amplification of genomic DNA:

1. VetPCR™ BAB.BOVI PCR Pre-mixture
2. DNase/RNase free water
3. BAB.BOVI PCR Positive control
4. BAB.BOVI PCR Positive control Pre-Mixture
5. Brig™ Molecular Weight marker
6. Mineral Oil Solution.

QUALITY CONTROL - This batch meets the following release criteria: Each components of the system has been tested and pased its pH, conductivity or RNase activity measurements acording to individual quality control specifications.

Bioingentech Genomic DNA Purification Kit is a registred trademark of Bioingentech Corporation.

Unless other stated, the performance of this product is guaranted for six months from the date of purchase if stored and handled properly.

Bioingentech Ltd.

Salas 350, Piso 2

Concepción

CHILE

Telephone (+56) 41-2790435

Website www.bioingentech.com

Info info@bioingentech.com

2004 – 2009 BioInGenTech corporation.
All rights reserved.