

## Product Information

**Kit Name:** Bioingentech Genomic DNA Purification Kit

**Storage Recomendations:** 22-25°C

This system is for preparation of high molecular weight genomic DNA from whole blood.

The purified genomic DNA is used for subsequent, i.e. amplification and restriction digest without further manipulation.

The Kit consist of the following components that are enough for purifying genomic DNA from 100x300µl whole blood.

1. Cell Lysis Solution
2. Nuclei lysis Solution
3. Proteinase K
4. RNase Solution
5. Wash Solution
6. Proteinase Precipitation Solution
7. DNA Rehydration Solution

**QUALITY CONTROL - This batch meets the following release criteria:**

Each components of the system has been tested and pased its pH, consuctivity or RNase activity measurements acording to individual quality control specifications.

**Bioingentech Genomic DNA Purification Kit** is a registred trademark of Bioingentech Corporation.

**Unless other stated, the performance of this product is guaranted for six months from the date of purchase if stored and handled properly.**

**PRODUCTS USE LIMITATIONS, WARRANTY, DISCLAIMER.**

Bioingentech manufactures products for a number of intendend uses. Please refer to the product label for the intendend use statements for specifics products. Bioingentech products contain chemicals which may be harmful if misused. Due care should be exercised with all Bioingentech products to prevent direct human contact.

Each Bioingentech products is shipped with documentation stating specifications and other technical information. Bioingentech products are warranted to meet or exceed the stated specifications. Bioingentech's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Bioingentech makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING WITHOUTH LIMITATION AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION OR ANY OTHER MATTER WITH RESPECT TO BIOINGENTECH PRODUCTS. In no event shall Bioingentech be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whethever based upon warranty, contract, tort (including negligence), or strict liability arising in connection with the sale or the failure of Bioingentech products to perform in accordance whit the stated specifications.

Bioingentech Ltd.

Salas 350, Piso 2

Concepción

CHILE

Telephone (+56) 41-2790435

Website [www.bioingentech.com](http://www.bioingentech.com)

Info [info@bioingentech.com](mailto:info@bioingentech.com)

2004 – 2009 BioInGenTech corporation.  
All rights reserved.